

1. One of the primary benefits of the MAXAIR CAPR Systems is:
 - a. Provides even airflow across the face
 - b. Does not need a battery belt
 - c. Provides greater freedom of movement
2. MAXAIR CAPRs can be used for protection from airborne particles. Examples may include:
 - a. Chemical and Biological
 - b. Chemical only
 - c. CBRN – Chemical, Biological, Radiological and Nuclear
 - d. Bacterial and Viral
3. Stethoscopes are able to be used with the MAXAIR CAPR cuffs because:
 - a. The cuff material is very thin and you can clearly hear through this material
 - b. The cuff of the DLC fits to the face in front of the ears
 - c. The airflow exits at the bottom of the cuff and is very quiet
4. MAXAIR CAPR helmets have 5 LED indicators to indicate:
 - a. 1 hour of battery life left and/or a damaged filter
 - b. Low airflow and/or amount of battery charge remaining
 - c. Battery is not functional and/or the filter needs immediate changing
5. The helmet must be worn with the front headband approximately ½” above the eyebrows because:
 - a. This allows for proper airflow, a wide field of vision, and easy visibility of the LED indicators in your peripheral vision
 - b. It allows for maximum air to be exhausted below the cuff or shroud
 - c. It fits better and allows less contaminant to enter the helmet
6. Adjustment of the MAXAIR CAPR helmet is accomplished by:
 - a. Adjusting the center band at the top of the helmet
 - b. There is no adjustment for the helmet as one size fits all
 - c. Adjusting the ratchet knob at the back of the helmet for circumference, and the helmet side headband adjustments for positioning on the head.
7. The airflow setting in the MAXAIR CAPR helmets from the factory is set at 6 CFM (cubic feet/minute)
8. The airflow adjustment switch allows you to change the airflow setting from:
 - a. 7 – 10 CFM
 - b. 6 – 9 CFM
 - c. 6.5 – 9.5 CFM
9. The reason there is no fit testing required with MAXAIR CAPRs is:
 - a. The motor runs more efficiently than any other PAPR on the market
 - b. There is greater area for airflow
 - c. They are positive pressure devices

10. The HE filter used with MAXAIR CAPRs provides what level of filtration efficiency:
 - a. 99.97% efficiency
 - b. Is rated with the same efficiency as that of an N95
 - c. 95% efficiency
11. Assembly of the HE filter cartridge onto the helmet is from:
 - a. Back snap then side snaps
 - b. Side snaps then back snap
 - c. It does not make any difference how the filter is assembled onto the helmet
12. The Lithium-Ion battery that is used with MAXAIR Systems for Emergency Preparedness typically provides for (hours per full charge):
 - a. 16-20+ hours
 - b. 8 hours
 - c. 4 hours and must be changed out with a new battery after lunch
13. The Lithium-Ion battery that is standard with MAXAIR CAPR for routine Infection Prevention applications typically provides for (hours per full charge):
 - a. 12 hours
 - b. 10 hours
 - c. 8 –10+ hours
14. The Lithium-Ion battery has no **MEMORY** which allows you to place the battery on the battery charger after each use, regardless of how long it was in actual use.
15. The Lithium-Ion battery should not be left on the charger without use for longer than:
 - a. it takes to become fully charged (<5 hours for small battery; <10 hours for large battery)
 - b. 1 week
 - c. Indefinitely
16. MAXAIR CAPR helmets and filter cover caps can be disinfected with the following:
 - a. Quaternary disinfectant wipes
 - b. Soap and water
 - c. Bleach diluted with water
 - d. Quaternary Ammonia, Bleach or Alcohol based disinfectant wipes